## **AT51005**

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**USER MANUAL ROLLATOR**

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**MANUEL D'UTILISATION ROLLATOR**

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**MANUAL DEL USUARIO ROLLATOR**

**MANUALE D'USO ROLLATOR**

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IU\_AT51005

### WSTĘP

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Rys. 3 Blokada

## INTRODUCTION

Thank you for purchasing our product. Before use, please read this user manual carefully. Incorrect use may lead to, injuries or accidents. Therefore, in order to use the rollator safely and successfully, please read this manual carefully.

**NOTE:** Inspect all parts for shipping damage. If any such damage is found, DO NOT use the product. More information you can get from the seller.

## GENERAL WARNINGS

Do NOT assemble and use the product until you have read and understood this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional, distributor, or service technician before attempting to assemble the product to avoid personal injury or damage to the product. Each user should ask his doctor or therapist to determine the appropriate variant of regulation and method of use.

The rollator is not intended to be moved by the user who is sitting on it.

A physician / therapist should be present when adjusting the height of the rollator to provide support for the user and to properly tighten the screws.

When the rollator is unfolded, always make sure that the tube under the seat is locked and properly attached to the frame. Make sure the height adjustment knobs are locked so that the main support structure will not move under pressure. Failure to follow the above recommendations may lead to injuries. Before use, make sure that ALL the adjustment elements of the handles and the seat are secure and that the wheels and moving parts of the product are in good condition.

During use, ALL the wheels must be in contact with the ground ALL THE TIME to ensure the stability of the rollator.

Failure to do so may cause the rollator to tip over and injure the user.

**CAUTION:** When using and operating the rollator, and when folding it and adjusting the mechanisms, there may be a risk of blocking and / or pinching of body parts of the user / accompanying person in openings / gaps between the components. These steps must be performed with extreme caution. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

### CAUTION:

When a product-related „serious incident“

occurs, which directly or indirectly led, could or may have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person, or
- c) a serious threat to public health

this „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

### CAUTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional.

**ATTENTION:** The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing or resulting from failure to comply with the recommendations contained in this manual.

**ATTENTION:** it is forbidden to use the product other than in accordance with its intended purpose

**NOTE:** If the product is used improperly, there may be a risk of tipping over. Please follow the using / moving instructions.

## WARNINGS REGARDING FOLDING AND ASSEMBLY

After assembling the rollator, make sure that the seat is locked and both front wheels are secured with the nuts. Also make sure that both handles are at the same height and do not move. BEFORE use, always check that the rollator and its components are securely locked in place.

**INTENDED USE:** The product is intended to support the movement of people with reduced walking ability. The rollator is designed for one person.

## CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

## PRODUCT FEATURES:

- Tilttable footrests, seat with backrest.
- Comes with a plastic shopping bag with a shoulder strap or attachment to the front of the rollator.
- Crutch / cane holder.

- Full wheels made of a material imitating rubber (40 mm wide), front swivel wheels, rear wheels with brakes.
- Brakes with the possibility of blocking.

**SCOPE OF DELIVERY:** ROLLATOR, BAG, INSTRUCTIONS

**UNFOLDING**

Remove front wheel hole covers, install wheels and tighten with nuts. Then mount the rear wheels, make sure the copper pins are in the correct holes. Then position the rollator and lock the tube under the seat. put the bag in it and then press the seat. Place the handles at the appropriate height (6 degrees of height adjustment - to adjust the height properly, it is recommended to consult a healthcare professional, tighten the knobs. Then mount the backrest so that the copper pins are on the outside. end to mount the footrests.

**FOLDING**

All disassembly operations should be performed in the reverse order to that described in the section „DISASSEMBLY“.

**APPLICATION OF HAND BRAKES**

Pull the brake towards the handgrips. Do one of the following:

1. STILL - hold the brake handle up (Fig. 2)
2. MOVING - release the brake handle (Fig.1)
3. TO LOCK - Depress the brake down.

**ATTENTION:** it is forbidden to use the product other than in accordance with its intended purpose

**STORAGE:**

The product should be stored in a dry place, not exposed to direct sunlight.

**MAINTENANCE**

1. Please ensure that all attachments are always secured.
2. Check all mounting bolts, tighten if necessary.
3. Replace any damaged parts immediately.
4. The product should be periodically cleaned with water and cleaning agents. The seat should be cleaned using a mild detergent.

5. When not in use, the product should be stored in a dry room. Avoid prolonged exposure to direct sunlight.

**WHEEL MAINTENANCE**

1. Wheels should be cleaned with warm water and a mild detergent. Dry with a dry cloth.
2. Lubricate periodically with 3-in-1 grease if the wheels begin to squeak.

**TECHNICAL PARAMETERS**

Foldable, four-wheeled aluminum rollator.  
 Wheel Size: 8”  
 Overall height: from 78 cm to 98 cm  
 Overall length: 81cm  
 Width in the upper part: 55 cm  
 Overall width: 60 cm  
 Seat width: 43 cm  
 Height from the ground to the seat: 56 cm  
 Seat depth: 29 cm  
 Distance between the handles: 49 cm  
 Weight: 9.7 kg.  
 Maximum user weight: 136 kg

This sign indicates the maximum weight of the

**TARGET PATIENT GROUP**

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

**DISPOSAL OF THE PRODUCT**

The product must be disposed of in accordance with the applicable disposal conditions.


Fig.1 Free movement

Fig.2 Braking

Fig.3 Lock

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WARRANTY CARD

Modell: .....  
Serial number: .....

Point of sale stamp  
and a legible signature of the Seller: .....

1. ANTAR MEDICAL Sp. z o.o. 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.  
03-068 Warsaw, ul. Zawiańska 43  
Tel. 22 518 36 00, Fax 22 518 36 30  
www.antar.net

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