AT51039

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USER MANUAL ROLLATOR

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NÁVOD NA POUŽITIE ROLÁTOR

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INTRODUCTION

Thank you for purchasing our product. Before using, please read this user manual carefully. Incorrect use may lead to health detriment, injury or accidents. Therefore, in order to use the product safely and successfully, please read these instructions carefully.

NOTE: Inspect all parts for damage caused in transport. If any such damage is found, DO NOT use the product. More information from the manufacturer.

RECOMMENDATIONS SAFETY

To be safe when using the walker, please follow these guidelines:

GENERAL WARNING

Do NOT assemble and use the product until you have read and understood this manual. If you do not understand any warnings, cautions, or instructions, contact your healthcare professional, distributor, or technician before attempting to assemble the product to avoid personal injury or damage to the product. Each user should ask his doctor or therapist to determine the appropriate variant of regulation and method of use. The walker is not intended to be guided by the user while sitting on it. A physician / therapist should be present when adjusting the height of a walker to provide support for the user and to properly adjust the brakes and tighten the screws. When the walker is unfolded, always make sure that the seat tube is locked and properly attached to the frame and that the brake cable is inside the side tube and is securely attached to it. Check the brake cable after each height adjustment to avoid loosening the cable. Make sure the knobs are locked so that the main support structure will not move as a result of pressure. The quick release mechanism is installed in place and will remain in the release position even if you pull forcefully in the opposite direction. The copper pins on the adjustment tube pop out from outside. Failure to follow the above recommendations may lead to injuries. Before use, make sure that ALL the adjustment elements of the handles and the seat are secure and that the wheels and moving parts of the product are in good condition. During use, ALL the wheels must be in contact with the ground AT ALL TIMES to keep the walker stable. BEFORE the user is seated on the seat, the brakes MUST be in the locked position. When the walker is used in a stationary position, its handbrakes MUST be in the locked position.

ATTENTION: The manufacturer is not responsible for damage caused by neglect of maintenance, inadequate servicing or resulting from failure to comply with the recommendations contained in this manual.

ATTENTION: it is forbidden to use the product other than in accordance with its intended purpose!

NOTE: In the event of a product-related „serious incident“ that led, directly or indirectly, could or could have led to any of the following events:

- a) death of a patient, user or other person or
- b) temporary or permanent deterioration in the health of a patient, user or other person, or

c) a serious threat to public health

this „serious incident“ must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE: In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of a medical device, please consult a healthcare professional.

NOTE: There may be a risk of tipping over if the product is used incorrectly. Please follow the sitting / standing up / moving instructions. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

CAUTION: When using and handling the product, as well as when assembling and adjusting the mechanisms, there may be a risk of entrapment and / or pinching of body parts of the user / accompanying person in openings / gaps between the elements. These steps should be performed with extreme caution. After completing the adjustment, stabilize the position by tightening the nuts / bolts carefully.

DESTINY

A rollator is intended for people with limited mobility who still have the ability, strength and balance to leave their home with the help of a rollator.

CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent the safe handling of the product.

WARNINGS ASSEMBLY AND ASSEMBLY


After disassembling or assembling the walker, make sure that all components are properly assembled.

INCLUDED

User manual, rollator, bag, cane holder.

ASSEMBLY

1. Raise the top.
2. Unfold horizontally, press the seat from above until you hear a click.
3. Install the cane / umbrella holder.
4. Assemble the bag.

FOLDING

Reverse the order of the ASSEMBLY section

APPLICATION OF HAND BRAKES

1. Pull the brake towards the handgrips.
2. Do one of the following:
 - A. STILL - hold the brake handle up:
 - B. MOVING - Release the brake handle


Fig. 5 (unlocked state)


Fig. 6 (locked state)

HANDBRAKE LOCKING

1. Push down on the bottom of the brake handle in compliance with the drawing until you hear a click (Fig. 6)
2. For more information, see UNLOCKING THE HANDBRAKE in this manual.

UNLOCKING THE HANDBRAKE

1. Pull the brake handles upwards until they are released as shown in fig. 5.
2. Release the brake handle.
3. Make sure the brake cable is placed along the side-line and the unused part is coiled in front of the walker.

TECHNICAL PARAMETERS

Overall width: 64 cm
 Overall length: 82 cm
 Overall height: 92.5 cm
 Seat width: 45 cm
 Seat depth: 23.5 cm
 Seat height: 55 cm
 Rear wheels: 8-inch PVC wheels
 Front wheels: 9-inch PVC wheels
 Weight: 8.1kg
 Package size: 86*26*37 cm
 Adjustable handle height: 79-91 cm
 Maximum payload: 136 kg

THIS SIGN MEANS A MAXIMUM
WEIGHT OF THE USER!

CLEANING

Rollator can be cleaned with a damp cloth or sponge with warm soapy water. Domestic anti-bacterial surface cleaners may be used if required. Care should be taken not to leave any cleaning residues on the handgrips or seat (after cleaning) which may make them slippery when used. Do not use scouring pads as this may damage surface coatings.

MAINTENANCE

1. Please ensure that all attachments are always secured.
2. Check the bolt securing the seat to the side elements. Tighten if needed.
3. Replace any damaged or damaged parts immediately.

WHEEL MAINTENANCE

1. Wheels should be cleaned with warm water and a mild detergent. Pat dry with a dry cloth.
2. Periodically lubricate with 3-in-1 grease if the wheels begin to squeak.

STORAGE

The product should be stored in a dry and cool room. Do not expose to direct sunlight.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

DISPOSAL OF THE PRODUCT

The product must be disposed of in accordance with the applicable disposal conditions.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłńska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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