## **AT51001**

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**WALKER USER MANUAL**

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**MANUALE D'USO DEL CAMMINATORE**

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**MD**

**CE**

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## INTENDED USE

The folding walker is designed for people with problems with verticalization, with diseases of the locomotor organs or people with problems with maintaining balance while walking. The walker provides the user with support, reducing the risk of tipping over. The product can be used when patient learns to walk or after injuries.

## CHARACTERISTICS

- made of reinforced aluminium
- equipped with a set of anti-slip feet
- easy to fold and transport, folded with one button
- 8 levels of height adjustment
- ergonomic foam grips

## CONTRAINDICATIONS

physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

## WARNINGS:

1. Before using the walker, read this instruction manual and keep it for future reference.
2. Using the walker on a ground that does not provide reliable support, i.e. wet, icy, greasy, heavily sandy, is unacceptable because of the risk of falling.
3. Before each use, - perform a visual inspection of all structural elements of the product, including: - carefully check of the stability and condition of the load-bearing frame
4. The walker can only be used to support one user; it is forbidden to stand on the walker.
5. Do not exceed the maximum height level.
6. Do not use if user weight exceeds 110 kg.
7. Do not use the walker if any visible or tangible damage is found (e.g. broken frame) as this may result in an accident.
8. Do not lean out of the walker. This may cause the user to fall over.
9. Do not swing the walker. This may cause

the user to fall over.

10. Avoid contact with open fire

## CAUTION:

In the event of a product-related „serious incident“ that directly or indirectly led to, may have led to, or may lead to, any of the following:

- a) death of a patient, user or other person or
  - b) temporary or permanent deterioration of the patient's, user's or other person's health or
  - (c) a serious threat to public health
- this „serious incident“ should be reported to the manufacturer and to the competent authority of the Member State where the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

## CAUTION:

In the event of pain, allergic reactions or other disturbing, unclear for the user symptoms related to the use of the medical device, consult a healthcare professional.

**ATTENTION:** it is forbidden to use the product in a way other than for its intended purpose!

**ATTENTION:** Before use, check that all elements have been assembled correctly.

**ATTENTION:** During the use and handling of the product, as well as during its assembly and adjustment of mechanisms, there may be a risk of entrapment and/or squeezing of the user's/assistant's body parts in the openings/gaps between the elements. These steps should be performed with particular care. After adjustment, stabilize the position by carefully tightening the nuts/bolts

**NOTE:** The manufacturer is not responsible for damage caused by neglected

maintenance, improper servicing or failure to comply with the recommendations contained in this manual.

### TECHNICAL PARAMETERS

Height: 75 - 93 cm, 8 levels of adjustment, every 2.5 cm

Width between handles: 46 cm

Width: 58 cm

Folded width: 9 cm

Length: 46 cm

Packaging dimensions: 52 x 10 x 79 cm

Maximum user weight: 110 kg

Weight: 2.2 kg

Weight with carton: 3 kg

THIS MARK INDICATES THE  
MAXIMUM USER WEIGHT.

### UNFOLDING

To unfold the walker, pull the sides of the walker outside until you hear the characteristic click of the latch, and then set the legs to the appropriate height.

### FOLDING

To fold the walker, press the latch button and pull the sides of the walker inwards, and then fold the walker.

### USE

Stand by the walker, holding the handles of the walker with both hands. Moving with this type of equipment consists in alternately lifting the sides of the walker up, moving them forward and then moving forward user's legs.

### TARGET GROUP OF PATIENTS

The healthcare professional shall indicate the use of the device to adults and children on his/her own responsibility, taking into account the available variants/sizes/necessary functions/sizes and indications, bearing in mind the information provided by the manufacturer.

### PRODUCT MAINTENANCE

The product should be periodically cleaned with water and mild detergents. The product may be subjected to action disinfectants. Make sure that the disinfectant does not contain harsh chemicals such as strong acids and alkalis. Do not use rough cleaning objects or strong detergents.

### STORAGE

When not in use, the product should be stored in a dry and cool place in its original packaging. Avoid damp rooms and prolonged exposure to direct sunlight.

### METHOD OF DISPOSAL OF THE PRODUCT AFTER WITHDRAWAL FROM USE

After the device is decommissioned, the medical device can be disposed of as ordinary municipal waste, with the exception of electrical products - proceed in a manner appropriate for the disposal of electrical and electronic equipment.

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WARRANTY CARD

Model: .....  
 Serial number: .....

Point of sale stamp  
 and a legible signature of the Seller: .....

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.  
 03-068 Warsaw, ul. Zawiańska 43  
 Tel. 22 518 36 00, Fax 22 518 36 30  
 www.antar.net

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