AT52320, AT52321

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ENG

INTRODUCTION

Before attempting to operate your power wheelchair for the first time, please read and observe all recommendations, warnings, and cautions in this manual. If any of the information in this manual is not understood, or if additional assistance is required for assembly or operation, please contact your dealer. The safe use of the product depends on compliance with the recommendations, cautions and warnings contained in this manual. We shall not be liable for any damage and/or injury resulting from individual unsafe operation or failure to observe the recommendations, cautions and warnings contained in this manual. These instructions are intended to ensure maximum safety. Understanding the instructions is an essential requirement for the safe operation of the wheelchair.

SPECIFICATIONS

Overall size(L*W*H): 1125*630*1305 mm
Seat height: 515 mm
Seat width: 440 mm
Seat depth: 450 mm
Height of armrests: 200 mm
Backrest height: 800 mm
Battery weight: 4.2 kg
Weight of wheelchair including battery: 40.1 kg
Maximum speed: 6km/h
Stopping distance at 6km/h: ≤ 1500 mm
Max. turning radius: ≤ 1200 mm
Payload: 120 kg
Theoretical range on a single charge: 28 km
Static stability: ≥ 9°
Dynamic stability: ≥ 6°
Climbing capacity: 6°
Motor specification: 24V/250W*2
Battery specification: 24V/20Ah*1
Maximum output current of the control unit: 2A
Battery output current: 30A
Front wheels: outer diameter 190 mm
Inflatable rear wheels: outer diameter 405 mm

DESCRIPTION

The powered wheelchair consists mainly of several parts:

Seat frame: consists of a folding component that can be folded from the centre of the seat base for convenient storage or transport. Operating system: consists of a control unit, a battery compartment, two motors, a set of drive mechanisms and drive wheels.

The battery compartment is located under the seat, housing a 24 V, 20 AH battery and its overload protection. The battery charger socket is located at the rear of the control unit. This power wheelchair is equipped with two manual slow-speed levers located on the drive mechanism. Adjusting the levers can stop the wheelchair.

DESCRIPTION

1. handle for the carer
2. backrest cushion
3. hand brake
4. pneumatic rear wheel
5. lane front wheels
6. anti-slip armrest
7. detachable footrest
8. calf rest
9. footrest plate
10. removable headrest
11. Rear additional control device (model AT52320)
12. remote additional control (model AT52321)

ASSEMBLY

1. Remove the trolley from the packaging and lay it on the ground. Unfold the frame of the pushchair on the left and right sides and push the two tubes on either side of the seat base into the grooves on either side of the frame of the pushchair (see fig. 2).
2. Lift up the handles and lock the folding device (see fig. 3).
3. Remove the screws and gaskets from the stem of the control unit, then insert the control unit and fit into the plastic gaskets, insert the screw and tighten with a screwdriver (see fig.5)
4. Then fix the wire harness to the frame with plastic cable ties (see fig.6).

Fig. 2

Fig. 3

Fig. 5

Fig. 4

Fig. 6

Fig. 7

Fig. 8

Fig. 9

Fig. 10

Fig. 11

5. Connect the control unit wiring harness plug to the battery wiring harness socket (see fig. 7).
6. Fit the footrest attachment to the top tube of the left or right frame (see fig. 8). And then rotate the footrests forward so that they are automatically locked, then unfold them (see fig. 9).
7. Attach the leg support strap to the right and left upper tube, then use the nylon strap.
8. Adjust the height of the footrest, unscrew the lock nut on either side of the footrest, (see fig. 11), move up and down the telescopic tube, adjust to a convenient height, then insert the bolts and tighten the nut.

USE

Control device

The control device is a key piece of wheelchair equipment.

All electrical components for operating the wheelchair are contained in the control device.


The control device consists of the following parts:

1. Joystick
2. Battery level indicator
3. On/Off button
4. Speed indicator
5. Increase speed button
6. Button for decreasing speed

If the wheelchair moves in an unexpected direction, release the joystick immediately, which will bring the wheelchair to a stop.

7. Horn button
8. Charger socket

The control unit is usually located on one of the armrests and is connected to the battery along with the motors.

On/Off button

Controls the power supply to the electronics of the control unit, which then drives the motors, do not use the on/off button to stop the wheelchair unless there is an emergency, otherwise the life of the drive components may be reduced.

Joystick

The joystick is mainly used to control the speed and direction of the wheelchair. The further you move the joystick from the central position, the faster the trolley moves. When the joystick is released, it automatically returns to the central position and automatic braking takes place.

WARNING


If the wheelchair accidentally moves, please release the joystick immediately so that the wheelchair stops moving automatically, unless the joystick is faulty.

Horn button

When this button is pressed, the horn will be activated.

Speed increase/decrease buttons

When the power is switched on, the speed light will indicate the maximum speed at which the wheelchair is currently moving. The range of maximum speed indicated by the speed light number can be adjusted by users. Each time the speed up (or speed down) button is pressed, the speed is increased or decreased by one unit.

Brake lever

If not using a wheelchair, pull the brake lever backwards to apply pressure to the wheels, which will immobilise the wheelchair.

If using a wheelchair (whether automatic or manual), push the brake lever forward to ensure that the lever does not press on the wheels (see fig. 12).

NOTE

When there is no need to drive on an incline, the brake lever must be moved forward, otherwise the wheelchair will lose steering and could lead to personal injury.

Seat belt

For your own safety, always fasten the seat belt. Engage the belt buckle so that you hear a click.

Charging the battery

The separate battery charger is the main part of the wheelchair. With it you can charge the wheelchair quickly and easily.

WARNING

The wheelchair battery must be charged using the charger provided. Do not use chargers designed for cars.

Charge the battery using the separate charger supplied with the wheelchair:

- Ensure that the control unit is switched off and the wheelchair is not in slow mode.
- Plug the 3-pin metal charger plug into the 3-pin charger socket on the control unit (see fig. 13).
- Insert the input plug of the charger into a wall socket.
- The lighting of the red LED on the charger indicates that charging is in progress. When the green LED on the charger lights up, it means that the battery is fully charged.
- The operating time is approximately 8-12 after one charge.
- When charging is complete, first remove the charger input plug from the wall socket and then the output plug from the socket on the control unit. Place the charger with cables and plugs together in the bag on the back of the seat back.

Fig.13

Overload protection

If the motors are overloaded, the overload protection will disconnect the power supply to protect the motors and their electrical components. Expert assistance is required to restore the protective function, which involves unscrewing the screw in the battery box, replacing the fuse, re-covering the battery cover, and then tightening the screw.

Precautions

Ensure that the control unit is securely mounted and that the joystick in the central position is upright.

Sit securely, fasten your seat belt, lean against the backrest, and raise your head to avoid injury while driving.

When the power is switched on, first check that the maximum speed indicated by the speed indicator light on the control unit corresponds to the user. If it does not, adjust it using the up or down speed button on the control unit. For safety reasons, we suggest starting at a lower speed and then gradually increasing the speed.

If the wheelchair moves in an unexpected direction, release the joystick immediately, which will bring the wheelchair to a stop.

Note

Keep your feet within the footrests, do not extend them out of the footrests.

Do not extend your hands from the armrests while driving.

Make sure the wheelchair is not in slow motion mode before operating or leaving the wheelchair.

Hazards

Do not use the wheelchair in the following cases:

1. cases prohibited in this manual, such as exceeding the maximum gradient, obstacles that are too high, motorways, motor vehicle lanes
2. in certain locations or on surfaces, such as wet grassy slopes, where the wheelchair may slip.
3. when the control system or other important components need to be repaired.

Electromagnetic interference

The wheelchair may run in areas affected by electromagnetic interference from certain radio transmitters such as radios, wireless intercoms, mobile phones, radar, etc. In these cases, the handling of the trolley may be affected.

WARNING! Electronic equipment can be exposed to electromagnetic interference (EMI). Such interference can come from radio stations, television stations, mobile phones, and other radio transmitters. If the wheelchair exhibits abnormal operation due to electromagnetic interference, switch off the power supply and consult a service technician. The company accepts no legal liability for losses caused by non-compliance with this condition.

INDICATIONS

The product is designed to comfortably transport people with limited or no walking ability. The wheelchair is designed to transport one person.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

FOLDING AND STORING THE WHEELCHAIR

To transport and store the wheelchair, fold it as follows.

1. place the wheelchair on the ground. Disconnect the power supply.
2. turn the footrests upwards and rotate outwards by 90° (see fig. 15).
3. pivot the footrest lever outwards, rotate it 90° and then pull it out (see fig. 16).
4. clamp the two folding levers on either side of the handle and lower the handle down (see fig. 17).
5. pull the seat base upwards by its centre and then fold the wheelchair (see fig. 18).

fig. 15

fig. 16

fig. 17



fig. 18

ADDITIONAL CONTROLLERS

AT52320

1. Connection cable
2. LCD display
3. Power on/off
4. Operating modes
5. Pivoting lever for forward and reverse travel
6. Wheelchair sign (illuminated in working condition)
7. Gear indicator
8. Battery charge indicator
9. Vertical mounting holes, M4.0 10.
10. horizontal mounting holes, m4.0
- 11., 12. Bolt holes

REAR CONTROL UNIT OPERATING INSTRUCTIONS

1. First install the attachment piece and then the rear control unit on the right-hand handle of the wheelchair, using the allen key to do so.
2. Plug the XLR/USB end of the rear device connection cable into the charging port (XLR or USB) on the wheelchair's main control unit and secure the cable with a cable tie.

SWITCHING THE REAR CONTROL UNIT ON/OFF

1. Hold down the "Power" button for 2s until the rear control device is switched on.
2. Press the "MODE" button to select the appropriate gear as required (gear range: 1-5).

USE

After steps 1 and 2, the operator of the rear control unit should hold the wheelchair handles with both hands and use his right thumb to control the pivoting lever. While pushing the lever forward, the wheelchair moves forward and vice versa; when the lever returns to the central position naturally, the electric wheelchair does not move.

FAULT INDICATION ALARM AND REMEDIES

Most common problems:

- (a) The wheelchair remains stationary regardless of whether the lever is moved forwards or backwards.
- (b) The wheelchair continues to move even when the lever is returned to the central position.

If these phenomena occur, check the parameters of the pivoting lever:

1. Disconnect the connecting cable from the main control unit.
2. Hold down the "Mode" button, meanwhile connect the rear control unit to the main control unit; when the rear control unit enters the verification mode, release the "Mode" button, and push the lever forward and backward at least twice (Note: push the lever as far as possible). Once the lever has returned to the central position, press the "Mode" button again to save the verified data and complete the process.

AT52321


1. Switch on/off
2. Horn
3. Release button
4. Control button
5. Acceleration button
6. RGP status

PAIRING THE CONTROLLER

Place two AAA batteries in the controller, then press the power button on the main.

When the RGP status indicates a blue light that flashes every 2s, press the horn and deceleration button.

Until the blue light is on continuously. The pairing has been correctly completed.

CLEANING AND MAINTENANCE

Care and maintenance after use

Switch off the power supply (preferably disconnect all cable connections). Do not allow children or unsuspecting persons to touch the wheelchair. Store the wheelchair in a cool and dry place to prevent it from reforming and to preserve its properties for longer. Use a clean and soft cloth to clean the frame and then allow to dry.

Routine checks

Before every ride, you should carry out an inspection of your trolley to keep it always in the best possible condition. In addition, inspect your trolley weekly, monthly and every six months according to the items in Table 2.

Troubleshooting

This trolley is equipped with an automatic alarm function that alerts you to malfunctions. If the trolley is not fully operational, the LEDs on the control unit will light up and a buzzer will sound to detect the fault according to the items in Table 3. If problems persist after troubleshooting, contact an authorised service centre.

Table 2 Basic troubleshooting

Number alert	LED display	Explanation
1	1 LED flashing	Low battery capacity
2	2 LED flashes	Left motor fault
3	3 LED flashes	Left brake fault
4	4 LED flashes	Failure of right motor
5	5 LED flashes	Failure of right hand brake
6	6 LED flashes	Control unit overload
7	7 LED flashes	Failure of the joystick
8	8 LED flashes	Failure of the control unit itself
9	9 LED flashes	Control unit failure

Table 3 Routine checks

Frequency	Always	Weekly	Monthly	Every 6 months
All parts			X	
Turning, driving, adjusting, and dismantling, etc.		X		
Brake	X			
Wiring harnesses and cables		X		
Battery	X			
Agility of the front wheels		X		
Pneumatic control box for the rear wheel		X		
Tire wear			X	
Wheel damage	X			
Handle, backrest, foot	X			
Engines				X
Control panel		X		
Cleanliness	X			

LIFETIME

We carried out a fatigue test covering the control unit, battery, motors, drive and frame of our trolley. The lifespan of the control device, drive and frame is on average 10 years, while the motor and battery are spare parts and therefore have a lifespan of up to 7 years.

The lifetime of a wheelchair in a clinical setting depends not only on the structural strength of the product itself, but also on the environment and the habits of the users of the product, whether periodic inspections, maintenance etc. are carried out. Therefore, the key to guaranteeing the longest possible service life is strict adherence by the user to this manual.

WARRANTY

Five year guarantee on the frame of the trolley.

Within one year from the date of purchase, for the following parts, we provide a free maintenance and replacement service to the original customer after the dealer has checked and found defects in materials and workmanship.

- Electric or lever control system
- Motor/drive system
- Shaft bearing and bushing

Six month battery warranty.

Please note that the warranty service is carried out by the dealer in conjunction with our after-sales department.

Outside warranty

- ABS plastic cover and rubber washer
- Tyres
- Decorative parts
- Damage due to abuse, mishandling, accidents and negligence
- Business use or other abnormal use

Please note that as part of care and maintenance, if you need to replace the battery, joystick, motor/drive, charger, and other service parts, please contact our company directly to replace the parts with new, original parts.

NOTICE

In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to or may lead to any of the following:

- (a) death of a patient, user, or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user, or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident.

In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

It is forbidden to use the product other than for its intended purpose.

Before use, check that all components have been assembled correctly.

The manufacturer will not be held responsible for damage caused by negligent maintenance, inadequate servicing or as a result of not following the instructions in this manual.

There may be a risk of tipping over if the product is used incorrectly. Please observe the instructions for getting on/off/rolling.

During use and handling of the product and when folding and adjusting the mechanisms, there may be a danger that parts of the user's/companion's body become trapped and/or pinched in the openings/gaps between components. Perform these operations with particular care. After adjustment, stabilise the position by tightening the nuts/bolts firmly.

This symbol indicates the maximum weight of the user!

HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE

After the device has been taken out of service, the medical device can be disposed of as normal household waste with the exception of electrical products - follow the procedures for the disposal of electrical and electronic equipment.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTA Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer.
2. During the warranty period ANTA undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTA reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTA (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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