AT52304

WÓZEK INWALIDZKI, ELEKTRYCZNY INSTRUKCJA OBSŁUGI
ELECTRIC WHEELCHAIR USER MANUAL
ELEKTOROLLSTUHL GEBRAUCHSANWEISUNG
UŽIVATELSKÁ PŘÍRUČKA K ELEKTRICKÉMU INVALIDNÍMU VOZÍKU
NÁVOD NA POUŽITIE ELEKTRICKÝ INVALIDNÝ VOZÍK
MANUEL D'UTILISATION DU FAUTEUIL ROULANT ÉLECTRIQUE
GEBRUIKERSHANDLEIDING ELEKTRISCHE ROLSTOEL
MANUAL DE USUARIO DE LA SILLA DE RUEDAS ELÉCTRICA
MANUALE D'USO DELLA SEDIA A ROTELLE ELETTRICA
ANVÄNDARHANDBOK FÖR ELEKTRISK RULLSTOL

IU_AT52304

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INTRODUCTION

Before using your power wheelchair, read and follow the information, warnings and cautions in this manual. If this manual contains information that you do not understand, or if you need further assistance with assembly or operation, please contact your dealer.

Following the instructions, warnings and cautions contained in this manual increases the safety of using the product. The manufacturer shall not be liable for any damage and/or injury resulting from improper use of the device or use not in accordance with the instructions and warnings contained herein.

APPLICATION

It is designed for people who cannot or do not want to move independently in a wheelchair using the strength of their own hands and for those who want to actively participate in everyday life.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling.

CAUTION! Failure to observe the warnings in the operating instructions may result in personal injury.

ATTENTION: Failure to follow the instructions in this manual may result in damage to the electric wheelchair.


ATTENTION: It is forbidden to use the product other than for its intended purpose!

DESIGN AND OPERATION

The electric wheelchair consists of the following main components (Figure 1): Seat frame: fitted with a mechanism that allows the seat to be folded along its axis, thus enabling convenient storage or transport. Control system: consisting of a controller, a battery compartment, two motors, a set of drive mechanisms and drive wheels.

The battery compartment is located under the seat. There is a 24V, 12AH battery with overload protection. The battery charger socket is located on the rear of the controller.

The drive mechanism is equipped with two manual levers. The levers allow the wheelchair to be stopped.

1.

TECHNICAL DATA

Overall size (Length*Width * Height): 1005*670*960 mm

Seat height: 520 mm

Seat width: 450 mm

Seat depth: 450 mm

Armrest height: 200 mm

Backrest height: 440 mm

Battery weight: 2.5 kg

Weight of wheelchair with battery: 36 kg

Maximum speed: 6 km/h

Braking distance at 6 km/h: ≤ 1500 mm

Max. turning radius ≤ 1200 mm

Payload: 120 kg

Approximate range on a single charge: 20 km

Static stability: $\leq 9^\circ$

Dynamic stability: $\leq 6^\circ$

Maximum angle of elevation: 6°

Motor specification: 24V/250W*2

Battery specification: 24V/12Ah*1

Maximum power supply to the controller: 50A

Battery supply current: 2A

Front wheel: non-inflated tyres, outside diameter 190 mm

Rear wheel: pneumatic tyres, outer diameter 406mm


THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

NOTE: It is forbidden to stand on the footrests as this can lead to a dangerous fall and injury.

MONTAGE

- 1 Remove the wheelchair from the packaging and place it on the ground. Pull the chair frame apart in opposite directions and press the two tubes on either side of the seat base down into the grooves on either side of the frame. (fig. 2)
2. Pull the handle upwards, blocking the folding of the wheelchair backrest. (fig. 3)
3. Connect the plug of the battery cable to the socket located in the frame under the seat. (fig.4)
4. Remove the screws from the controller base, then mount the controller to the bracket at the armrest and tighten the screws. (Figure 5).
5. Connect the controller wires marked R and L to the right and left motors respectively.
6. Then fix the wiring harness to the frame with plastic ties. (Fig. 6)
- 7 Fit the footrests to the left and right sides of the frame (fig. 8) Then turn the footrests over and unfold. (fig. 9)
8. The leg support strap should be placed over the right and left sleeves and then secured with tape.
9. Adjust the height of the leg rests, remove the lock nut on either side of the leg rest, remove the bolts (see fig. 11), adjust to the desired height, then insert the bolts and lock with the nut.

METHOD OF APPLICATION

CONTROLLER

The controller is a key component of the wheelchair. The controller contains all the components necessary for control (Figure 12).

The controller consists of the following parts:

- 1. handle
- 2. battery level indicator
- 3. on/off button
- 4. speed indicator light
- 5. acceleration button
- 6. speed reduction button
- 7. audible signal button
- 8. charger socket

12

The controller is usually located on one of the armrests and is powered by the batteries used to power the motors.

Controller switch-on button (see Fig. 12)

On/Off button

The On/Off button is responsible for controlling the power supply to the controller electronics and for controlling the power supply to the motor. Do not use the On/Off button to stop the wheelchair except in emergency situations. Using the On/Off button for this purpose may adversely affect the life of the product.

Control lever

The handle is used to control the speed and direction of the wheelchair. Speed is controlled by the degree of deflection of the handle relative to the centre axis. When the handle is released, the control returns to the centre position and the brake is automatically applied.


NOTE: In the event that the wheelchair moves on its own, it is essential to release the handle, causing the device to stop.

BEEP BUTTON

Pressing the button triggers an audible warning signal.

ACCELERATION BUTTON

When the power is switched on, the speed light will indicate the maximum speed of the wheelchair. The speed reached is displayed by LEDs and can be reduced by the user. Each time the speed up (or speed down) button is pressed, the speed indication displayed by the LEDs is increased (or decreased).

BRAKE LEVER

Move the brake lever backwards when starting the wheelchair so that the wheels are locked to prevent accidental movement. To move, move the brake lever forward so that the wheels are unlocked. (figure 13.)

13

NOTE: Be sure to use the brake only when necessary (e.g. driving on inclined surfaces), otherwise the brake may cause you to lose control of the wheelchair and cause a health risk to the user.

SAFETY BELT

To ensure your safety, remember to wear your seatbelt every time you use the wheelchair. Engage the belt buckle, noting the distinct clicking sound (Figure 14)

14

BATTERY CHARGING.


An external charger is part of the wheelchair's equipment. The charger enables the battery to be charged quickly and easily, ensuring its trouble-free operation.

NOTE: Wheelchair batteries should only be charged using the charger provided. Do not use chargers for vehicle batteries. Before first use, charge the battery fully first, not exceeding an 8-hour charging cycle.

THE WAY THE BATTERIES ARE CHARGED:

- Ensure the control panel is off and the wheelchair wheels are locked.
- Plug the charger plug into the socket located on the control panel (Figure 15).
- Insert the charger plug into the wall socket.
- The start of charging is indicated by a red LED on the charger. The lighting of the red, yellow and 3 green LEDs indicates that the battery is fully charged.
- The duration of a full charging cycle is between 8 and 12 hours.
- When charging is complete, first disconnect the charger plug from the wall socket and then the other plug from the driver socket. When charging is complete, the charger with cables and plugs should be placed in the bag located at the back of the backrest.

In the event of overloading the motors, there is an overload protection system that cuts off the power supply to protect the motors and the electrical installation. Restoring operation once the protection system has been activated requires the intervention of a specialist - remove the screws located on the battery compartment, replace the fuse and then reapply the battery compartment cover.

SAFETY BELT

To ensure your safety, remember to wear your seatbelt every time you use the wheelchair.

PRECAUTIONS

Ensure that the controller is correctly installed and that the handlebar is in the central position, pointing vertically upwards. The user should sit in the wheelchair, fasten the seat belt and rest the back against the backrest.

When the power is switched on, first check that the maximum speed indicated by the controller lamp is suitable for your needs. If it is necessary to change the maximum speed, adjust it using the up or down speed button on the control panel. For safety reasons, it is recommended to start adjusting the maximum speed from lower values. In the event that the wheelchair starts to move accidentally, immediately release the handle, which will stop the engine.

Note: The user's feet must be on the footrests.

Note: Keep your forearms on the armrests when using the wheelchair.

Note: Make sure the wheels are locked before using or putting away the product.

DANGER

The use of a wheelchair is prohibited:

1. When expressly forbidden in this manual, i.e. driving on steep gradients, overcoming excessively high off-road obstacles, driving on motor vehicle traffic routes, etc.

2. in certain locations or on surfaces that present hazards (e.g. excessive ice, dampness).
3. where a fault is found in the control system or other key components.

ELECTROMAGNETIC INTERFERENCE

The wheelchair may move through areas affected by electromagnetic interference caused by certain radio transmitters, such as radios, wireless intercoms, mobile phones, radars, etc. However, electromagnetic interference may affect the operation of the wheelchair.

NOTE

Electromagnetic interference (EMI) may affect the operation of the unit. If the unit is found to be malfunctioning due to the influence of electromagnetic interference, contact service immediately. The manufacturer will not be held liable for any damage resulting from failure to comply with the above recommendation.

FOLDING AND STORING THE WHEELCHAIR.

Store the wheelchair in a dry and cool room away from direct sunlight.

For transport or storage, fold the wheelchair according to the instructions indicated below:

1. place the wheelchair on the ground. Switch off the power supply. Raise the foot rests and rotate 90° (fig. 17).
2. unfold the leg rest section and turn it 90° and pull it out completely (fig. 18).
3. clamp the wheelchair lock release levers located on both sides of the handles turn (fig.19).
4. then pull the seat base out of the centre and fold the wheelchair (fig. 20)


17

18

19

20

SCOPE OF SUPPLY

Wheelchair, controller, 2 footrests, calf stabilisation belt, instruction manual.

GUARANTEE

1.1 Serial number

For after-sales service and warranty purposes, the vehicle identification number must be entered.

MODEL: AT52304

SERIAL NUMBER:

1.2 Content of the warranty

1.2.1 Repairs are carried out free of charge if the problems are of a manufacturing nature or are related to a (a repair of this kind will hereafter be referred to as a warranty repair).

1.2.2 Warranty repairs can only be carried out by AN TAR service.

1.2.3 AN TAR is not responsible for transport costs related to repair and replacement.

1.2.4 Replaced parts are the property of AN TAR.

1.3 Guarantee period and parts covered by the guarantee

Electrical parts (motor, controller, charger): 12 months from the date of purchase

Wheelchair frame: 12 months from date of purchase

Batteries: 6 months from date of purchase

1.4 Parts not covered by the warranty

The items listed below are not covered by the warranty:

1.4.1 Wear and tear of components resulting from use and the passage of time (e.g. normal material ageing, deterioration of quality of paint and plastics, etc.)

1.4.2 Damage caused by lack of maintenance or improper handling and storage;

1.5 Problems resulting from the following causes are not covered by the warranty:

1.5.1 Improper storage;

1.5.2 Interference by changing or adding non-original components to the vehicle;

1.5.3 Damage caused by the use of parts other than AN TAR original parts or other parts not recommended by AN TAR;

1.5.4 Excessive or improper use of the vehicle, such as driving over kerbs and overloading the wheelchair;

1.5.5 Driving in areas where the vehicle should not normally be engaged or used;

1.5.6 Damage or corrosion of surfaces caused by environmental factors such as precipitation air (chemicals, tree sap, etc.), stones, hail, windstorms, floods, etc.;

1.5.7 Modification of any materials or components not authorised by AN TAR;

1.5.8 The wheelchair may not be pulled or pushed by another vehicle when the wheelchair is occupant or the rear wheels are not set in "free" running mode.

CARE AND MAINTENANCE

1. Care and maintenance after use

Switch off the power supply (it is recommended to disconnect all electrical connections). The wheelchair should not be operated by children or persons with cognitive impairments. The wheelchair should be stored at room temperature so as to reduce any potential deformation caused by temperatures, thus ensuring that the device functions efficiently for a long period of time. Clean the wheelchair with a soft cloth and detergent, making sure to dry the device every time. Rough objects and strong detergents must not be used for cleaning.

2. Checking the correct functioning of the wheelchair

Several checks should be carried out before each ride to maintain the efficiency of the wheelchair. In addition, periodic checks are also recommended - once a week, once a month and every 6 months, as shown in Table 2.

3 Troubleshooting

The wheelchair is equipped with an automatic fault notification system. Malfunctions are indicated by LEDs on the controller and an audible signal. Malfunction information should be read in accordance with Table 3. If rectification of the fault does not stop the alarm signal, contact an authorised service representative.


Table 2 Error messages

| ERROR NUMBER | LEDS | DIAGNOZE |
|--------------|--------|-------------------------------|
| 1 | 1 LED | Low battery charge |
| 2 | 2 LEDs | Left engine failure |
| 3 | 3 LEDs | Failure of the left brake |
| 4 | 4 LEDs | Failure of the right engine |
| 5 | 5 LEDs | Failure of the right brake |
| 6 | 6 LEDs | Controller overload |
| 7 | 7 LEDs | Failure of the control lever |
| 8 | 8 LEDs | Controller system failure |
| 9 | 9 LEDs | Failure of the control system |

Table 3 Unit checks

| Controlled element | Daily | Once a week | Every month | Every six months |
|--|-------|-------------|-------------|------------------|
| Each part | | | X | |
| Turning, guiding, driving, setting, folding, etc. | | X | | |
| Braking | X | | | |
| Connecting cables | | X | | |
| Batteries and charging | X | | | |
| Front wheels | | X | | |
| Rear wheels | | X | | |
| Tyre wear | | | X | |
| Tyre damage | X | | | |
| Degree of wear of push handles and backrest and seat | X | | | |
| Engines | | | | X |
| Controller | | X | | |
| Cleanliness | X | | | |

ELECTRICAL DIAGRAM

IU_AT52304
 III ΔT52304

NOTE:

For intensive use of the product, it is recommended to replace the batteries/batteries every 6 months. The guarantee for the batteries/batteries is 6 months.

NOTE: Check rear tyre pressure (2.4-3bar) before first use. The pressure should also be checked periodically.

NOTE:

In the event of a product-related 'serious incident' that directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) the death of a patient, user or other person; or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

this 'serious incident' should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTE:

If you experience pain, allergic reactions or other distressing, unclear symptoms related to the use of the medical device, consult a healthcare professional.

WARNING: During use and handling of the product and when assembling and adjusting the mechanisms, there may be a danger of trapping and/or pinching parts of the user's/companion's body in the openings/gaps between components. Perform these operations with particular care.

NOTE: There may be a risk of tipping over if the product is used incorrectly. Please follow the instructions for getting on/off/ moving. Once the adjustment is complete, stabilise the position by tightening the nuts/bolts carefully.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

HOW TO DISPOSE OF THE PRODUCT

Dispose of the electrical medical device in accordance with the obligatory rules. Handle in a manner appropriate to the disposal of electrical and electronic equipment.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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