AT53043

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ELBOW ORTHOSIS USER MANUAL

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EN

GENERAL INFORMATION

The elbow orthosis is made of high-quality fabric that combines the advantages of nylon, polyester and latex, providing excellent strength, flexibility and comfort, as well as a proper fit. The orthosis is made of elastic material, using a flat weave knit fabric.

INTENDED USE

The orthosis provides strong and effective support for weakened muscles and the elbow joint and prevents injury. Silicone pellets will provide additional compression to the elbow joint. The product is designed to provide maximum support while maintaining a comfortable fit. The orthosis is made of flexible material, using a flat-weave knit fabric. It is important to fit the correct size to ensure maximum benefit and avoid excessive compression. Adult supervision is required for use by children. Keep out of reach of small children.

INDICATIONS

For support/stabilization and compression of the joint in degenerative, post-traumatic or post-operative soft tissue injuries in the elbow joint area, such as in the case of

- Inflammation of the tendon sheath in the elbow (tennis elbow)
- Elbow weakness
- Preventive measure to relieve and protect against injury and strain.

CONTRAINDICATIONS

Physical or mental limitations that prevent safe handling.

FEATURES:


silicone pellets provide additional increased compression in injury-prone regions
Extremely lightweight, material that provides "breathability" under the orthosis
comfortable hem

PRODUCT CLEANING INSTRUCTIONS:


- The orthosis can only be washed by hand in lukewarm water (up to 40°C) with a mild detergent (such as gray soap).
- Do not iron.
- Do not dry wash.
- Do not use bleach.
- Protect from moisture.
- Do not expose to high temperatures.
- Do not tumble dry.
- Allow to dry unfolded, squeeze gently before drying.

USE

The orthosis is worn by pulling it onto the hand at the height of the elbow joint. The product should be placed precisely on the elbow joint, so that the silicone pellets are on both sides of the elbow, stabilizing the elbow joint and the dark strip is on the inside of the hand.

To choose the right size, it is recommended to measure the circumference of the hand 10 cm under the elbow, in the natural straight position of the hand.

SIZES AVAILABLE: M (20-24cm), L (24-27cm), XL (27-30cm), XXL (30-33cm).

NOTES:

In the event of a device-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

the above "serious incident" must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

WARNINGS: do not use the product if warming fluids or ointments have been applied, do not use in case of skin damage.

WARNINGS: it is forbidden to use the product other than in accordance with its intended use.

NOTES:

- Before putting on the orthosis, read the contents of this manual.
- Before each use, make sure that there are no seam abrasions on the orthosis.
- Use of the product contrary to the rules contained in the instructions may result in serious complications.
- The orthosis should be used as recommended by a doctor or physiotherapist.
- The product should not be used in case of open wounds, abrasions of the epidermis, as well as after previous application of warming gel or ointment.
- Wearing the orthosis at night is possible only in the case of such a recommendation by the attending physician.
- The duration of wearing the orthosis should be consulted with the doctor.
- If soreness or rash occurs, discontinue use immediately and consult a doctor.
- If the orthosis is used after an injury, consult a doctor if the condition does not improve within 72 hours.
- If the injury is severe or if you have any doubts, consult a doctor immediately.

SUPPLY SCOPE: orthosis, instruction manual.

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this device is intended (see the section on the purpose of the device of this manual). The device may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the device on its own or as indicated by a doctor/therapist/other specialist, the available sizes/necessary functions and variants of the device, indications and contraindications for use, and information provided by the manufacturer should be taken into account.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

METHOD OF DISPOSAL OF THE PRODUCT AFTER WITHDRAWAL FROM USE

The product can be disposed of as ordinary municipal waste.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiślańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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