INSTRUKCJA UŻYTKOWANIA PNEUMATYCZNY SYSTEM TERAPII KOMPRESYJNEJ
USER MANUAL PNEUMATIC COMPRESSION THERAPY SYSTEM
GEBRAUCHSANWEISUNG PNEUMATISCHES KOMPRESSIIONSTHERAPIESYSTEM
NÁVOD K POUŽITÍ PNEUMATICKÝ SYSTÉM KOMPRESNÍ TERAPIE
NÁVOD NA POUŽITIE PNEUMATICKÝ SYSTÉM KOMPRESNEJ TERAPIE
MANUEL D'UTILISATION SYSTÈME DE THÉRAPIE PAR COMPRESSION PNEUMATIQUE
GEbruiksaanwijzing PNEUMATISCH COMPRESSIETHERAPIESYSTEEM
MANUAL DE USUARIO SISTEMA DE TERAPIA DE COMPRESIÓN NEUMÁTICA
MANUALE D'USO SISTEMA DI TERAPIA COMPRESSIVA PNEUMATICA
BRUKSANVISNING PNEUMATISKT KOMPRESSIONSBEHANDLINGSSYSTEM

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The compression therapy system consists of an air pressure sensor, air pump and cuffs, working together as one unit. The air pump is connected to the cuffs by special hoses. The direction of compression massage is from the extremities to the center of the body by sequentially inflating and then deflating the air chambers. The pressure can be adjusted to avoid patient discomfort. The cuff operates under the influence of a sensor and microprocessor. The compression therapy system improves blood and lymph circulation, preventing venous thrombosis and relieving lymphedema.

PURPOSE

The product is indicated for use by patients at home under medical supervision to reduce edema and prevent limb thrombosis or to treat lymphedema, post-traumatic and sports edema, edema after immobilization or venous insufficiency. The system can be used in patients over the age of 18 at home or in the hospital.

CONTRAINDICATIONS

1. acute inflammatory skin diseases
2. arrhythmia
3. rubella
4. deep vein thrombosis
5. pulmonary edema
6. acute venous thrombosis
7. fluctuating blood pressure
8. patients with artificial pacemakers installed
9. arterial insufficiency

SIDE EFFECTS.

1. discomfort if outlet pressure is too high
2. redness, itching
3. feeling of discomfort and pain during the procedure with a pressure of 120 mmHg

SAFETY MANUAL

It is important to read all the warnings and cautions in this manual, as they are intended to protect you, prevent injury, and avoid situations that may cause damage to the device.

- Do not use this device if you have an implanted pacemaker, defibrillator or any other metal or electronic device. Using such a device may cause electric shock, burns, electrical interference or death.
- Not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.
- The devices are not intended for personal unsupervised use by patients who are emotionally disturbed, suffer from dementia or mental illness.
- Repair is not allowed. If you find that the device is malfunctioning, do not disassemble it unwillingly or by yourself.
- If you feel abnormal or experience pain when using the product, stop using it immediately.
- Before the first use, consult a doctor, and the treatment should not last longer than 20 minutes. Thereafter, gradually increase the treatment time as recommended by the doctor (prolonged use may cause side effects).
- Keep the device away from infants and small children.
- The device should not be placed in areas where it is difficult to disconnect the power supply.
- Use only 220-240V power supply. Disconnect the device properly from the mains after treatment.
- Incorrect connection to the power supply may cause damage to the device. Keep your hands dry when connecting or disconnecting the power supply. Wet hands may cause electric shock
- The power cable must not be subjected to excessive stress, and no objects must be placed on it.
- Do not bend or fold the power cord. Damage to the cable may cause electric shock.

NOTES:

In the event of a product-related "serious incident" that directly or indirectly led to, may have led to, or may lead to any of the following:

- (a) death of a patient, user or other person or
- (b) temporary or permanent deterioration of the health of the patient, user or other person or
- (c) a serious threat to public health

the above "serious incident" should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient resides. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

If you experience pain, allergic reactions, or other worrisome, unclear symptoms related to the use of the medical device, consult your healthcare professional

ATTENTION: it is forbidden to use the product other than for its intended purpose.

NOTE: Before use, check that all components have been assembled correctly.

ELECTROMAGNETIC COMPATIBILITY (EMC) INFORMATION

1 This product requires special precautions for electromagnetic compatibility and must be installed and put into service in accordance with the following electromagnetic compatibility information. The product may be affected by portable and mobile RF communication equipment.

(2) Do not use cell phones or other devices that emit electromagnetic fields near the product. This may cause the product to malfunction.

3 This product has been thoroughly tested and verified to ensure proper operation.

4 The device should not be stored side by side or used with other devices, and if side-by-side use is required, the device should be observed to verify proper operation in the configuration in which it is used.

SPECIFICATION

Number of chambers: 4 chambers

Device size: 240×200×120 mm

Weight (with accessories): 2.22 kg

Treatment time: 10~30 minutes

Pressure during treatment: 30~250mmHg

Cycle time: 85s ~ 133s, ± 1s

Treatment modes: 1 mode

Input voltage: AC 220-240V, 50-60Hz

Power consumption: 65VA

Noise level: ≤ 55dB

Protection against electric shock: Class II, application part Type BF

Operating mode: continuous operation

Waterproof: IP21

Lifespan of the product: 5 years

Operating environment: +10 °C - + 40 °C, 10% - 95% 86kPa-106kPa

Storage conditions: -40 °C - + 70 °C, 10% -100% RH, 50kPa-106kPa

LEG CUFF SIZE (XXXL)

Total cuff length: 110 x 70cm

Recommended for people with height: 170-185cm

Maximum leg circumference: 62cm

HAND CUFF SIZE (L)

Total cuff length: 90 x 60cm

Recommended for people with height: 170-185cm

Maximum arm circumference: 52cm

BELLY CUFF SIZE (M)

Total waist length: 92-120cm x 37cm

Recommended for people with height: 170-185cm

Circumference (at widest point) : 80-108 cm

Instructions for Use:

1. connect the power cord.

2. connect the tube to the cuff(s) and place it on the body.

- To avoid knee pain (for leg cuffs), do not bend the cuffs during use

- It is not recommended to wear cuffs on bare skin

- Relax your muscles before application

- Close the zipper (if the cuff is equipped with it)

- Adjust the air outlets of the cuffs upwards (leg/hand cuff) and then connect the hose connector to the air outlet of the device

3. Turn on the device, then set the treatment time and pressure level. The standby state is displayed as a green light. The operating status is displayed as a blue light).

INSTRUCTION MANUAL

1 Place the device on a flat surface.

2. installation only requires connecting the connector to the device and cuff, and then connecting the mains power supply.

Note: When you receive the unit, open the filter cover and remove the security screw before use. The screw is used ONLY to secure and protect the unit from damage during transportation.

1 Plug the main plug of the power cord into a 220~240V outlet.

2 Check if the hose is well connected to the cuff(s): if yes, go to step 3; if no, insert the end of the hose into the outlet cuff(s).

3. connect the main end of the connection hose to the outlet of the main unit.

4 Check and make sure all connections to the main unit and the cuff are tight.

5. close the zipper at the cuffs.

Note: The connecting hose must not be worn around the neck to avoid strangulation or injury.

(6) Press the button to turn on the device.

7. set the pressure by pressing the button

Note: When using it for the first time, it is better to start with a pressure that is not too high, such as 120 or 150 mmHg, and then adjust the pressure to a level that is comfortable for the user.

8 Set the time by pressing the button.

Note: The time range is from 10 to 60 minutes. The default value is 10 minutes. When the time runs out, the device will turn off automatically.

Note: Do not open the cuff(s) during treatment. This may damage the cuff(s). Always turn off the device before removing the cuff(s).

9. after the cycle is completed, simply disconnect the tube from the device, the wires can be connected to the cuffs.


PRODUCT DESCRIPTION

1. The main unit - the product mainly consists of the device, cuff and connecting hoses.

NOTE: Before the first use, the user should open the filter cover and unscrew the screw located at the bottom of the device.

2 Control panel

- Time setting: dial to set the treatment time from 0-30 minutes
- Continuous operation: turn the knob to "Continuous operation" - the device runs all the time.
- Pressure setting: pressure range 30-250 mmHg, press the „plus“ button once to increase the pressure and press the „minus“ button once to decrease the pressure
- ON/OFF switch: press, the device is in standby mode; press again to turn off the device.

CLEANING AND MAINTENANCE

The following daily maintenance requirements must be fully adhered to in order to ensure the proper operation of the product and guarantee its long-term performance and safety.

CLEANING

Do not clean the device during treatment, make sure the device is unplugged before cleaning.

Clean the device with a soft, slightly damp cloth. For heavier dirt, you can also add a mild cleaner.

Make sure that water does not get into the device. Do not expose the device to moisture. Do not, under any circumstances, hold the device under running water or immerse it in water or other liquids. For hygienic reasons, each user should use his own set of cuffs. Do not use solvents or strong detergents to clean the cuffs. Do not wash the cuffs or clean them under running water. Gently clean the surface of the cuffs with a damp cloth and mild detergent.


CONSERVATION

1. if you have problems with the device, contact your dealer.
2. The manufacturer is not responsible for the consequences of repair work done on its own.
3. The user must not attempt to repair the device or its accessories. Contact your dealer for repair.
4. It is not allowed to open the device by unauthorized persons, this will void the warranty.
5. Each product in production has undergone systematic validation. Performance is stable and does not require calibration or validation.
6. If your product cannot achieve the expected performance, and the basic function has changed during normal use, contact your dealer.


PROBLEM SOLVING

Frequent faults	Check and confirm
If the device cannot start (if the plug is well connected to the power supply)	- whether the plug is well connected to the power supply - when the device is on - if the power supply is normal (220-240 V) - if the device shuts down
There are leaks and other strange noises	- if the connection hose and spigot are damaged - whether the plug is well connected to the device - if the connecting hose is crushed or kinked


LEG CUFF ARM CUFF

ARM CUFF

MANKEIT - BELLY BELT

TARGET PATIENT GROUP

Persons struggling with diseases, dysfunctions or injuries, for the treatment/rehabilitation or compensation of which this product is intended (see the section on the purpose of the product of this manual). The product may be purchased by the user independently or on the recommendation of a physician, therapist or other specialist. Whether purchasing the product on your own or as indicated by a doctor/therapist/other specialist, you should consider the available sizes/necessary functions/sizes/variants of the device, indications and contraindications for use, as well as the information provided by the manufacturer.

COMPONENTS

Device, hose set, cuffs (two leg cuffs are included), other cuffs can be purchased additionally.

STORAGE

Do not lead to sharp kinks in the connecting hoses or cuffs. After use, put the cuffs in the original package. Do not expose the device to direct sunlight, protect it from dirt and moisture. Store the device in a cool, well-ventilated place. Never place heavy objects on the device.

METHOD OF DISPOSAL OF THE PRODUCT

The medical device is an electrical product, please proceed in a manner appropriate for the disposal of electrical and electronic equipment.

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
BEINMANSCHETTE

Technical drawing of a mechanical part, likely a propeller hub, showing a side view and a cross-section. The side view shows a cylindrical base with a diameter of 19 cm, a total length of 90 cm, and a series of internal features with diameters of 23 cm, 24 cm, 25 cm, and 26 cm. The cross-section shows a 30 cm diameter and a 23 cm diameter feature.

Technical drawing of a rectangular plate with dimensions and internal structure. The plate has a total width of 120CM and a total height of 37CM. The width is divided into four equal sections of 30CM each, labeled 2, 3, 4, and 5. The height is divided into two main sections: a top section of 10CM and a bottom section of 27CM. The top section contains four horizontal rows of small circles, with a vertical dimension of 5CM between the top and bottom rows. The bottom section contains four horizontal rows of small circles, with a vertical dimension of 4.5CM between the top and bottom rows. A dashed line indicates the top edge of the plate.

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
MANGUITO DE PIERNA

Technical drawing of a mechanical part, likely a propeller hub, showing a side view and a cross-section. The side view shows a cylindrical base with a diameter of 19 cm, a total length of 90 cm, and a series of internal features with diameters of 23 cm, 24 cm, 25 cm, and 26 cm. The cross-section shows a 30 cm diameter and a 20 cm height.

Technical drawing of a rectangular box. The overall dimensions are 120CM in width and 37CM in height. The box is divided into four vertical sections of 10CM, 9CM, 9CM, and 9CM width. The rightmost section is labeled '1'. The top and bottom edges of the rightmost section are labeled '5CM'. The internal structure shows horizontal lines and circles, with labels '4.5CM' and '5CM' indicating specific dimensions. A dashed line indicates a fold or assembly line.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

FOSHAN HONGFENG CO., LTD
Add. No. 4-2 Legiang Road,
Leping Sanshui Foshan,
528100, Guangdong, Chins

Authorized Representative:
MedNet EC-REP GmbH
Borkstrasse 10 48163
Muenster Germany

Importer:
ANTAR MEDICAL Sp. z o.o.
03-068 Warszawa, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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